

Padycare[®] washable Incontinence products

Absorbing volume

	Legal requirements	Achived volume
Füssen	150 ml	266 ml
Prien	150 ml	469 ml
Starnberg	600 ml	1609 ml
Slip <u>Rügen und Sylt</u> Incl. absorbing core		
Size S	500 ml	739,0 ml
Size M-L	750 ml	1025,0 ml
Size XL-XXL	1000 ml	1526,5 ml
Slip Fehmarn Incl. absorbing core		
Size S	500 ml	841,0 ml
Size M-L	750 ml	962,8 ml
Size XL-XXL	1000 ml	1248,0 ml
Diaper Lindau		
Size S	600 ml	1492 ml
Size M	600 ml	1753 ml
Size L	600 ml	2037 ml
Size XL	600 ml	2425 ml
Size XXL	600 ml	3179 ml
Diaper Hof		
Size S	1000 ml	2196 ml
Size M	1000 ml	2733 ml
Size L	1000 ml	3073 ml
Size XL	1000 ml	3696 ml
Size XXL	1000 ml	4841 ml
Bed pads		
without wings	371,3 ml	1773 ml
with wings	477,4 ml	2242 ml